

PART LIST				
ITEM	QTY.	PART NUMBER (IMPIERIAL)	PART NUMBER (METRIC)	DESCRIPTION
1	1	-	PS-S	Pole stand
2	1	-	ZAS-FAC-PST	Focus unit
3	1	-	TSD-602C	XY Axis Steel Extended Contact Translation Stages
4	1	-	OUCI-CFH	C Mount Filter Unit
5	1	-	OUCI-BKT	OUCI Mounting Bracket
6	1	-	OUCI-3(NIR)_Cmount-sepa	Observation unit with coaxial illumination (NIR) Wavelength:460-1500nm
7	1	-	PAL-NIR	Objective Lenses
8	1	-	FHL-102	Wavelength:460-1500nm
9	1	-	MSL-1000F	Light Guides
10	1	-	STC-LBS132U3V-SWIR	SWIR camera Wavelength:460-1700nm

