



PART LIST				
ITEM	QTY.	PART NUMBER (IMPIERIAL)	PART NUMBER (METRIC)	DESCRIPTION
1	1	-	TAM-602SDM	60 mm XY Crossed Roller Stage with Diff Mics, +/- 6.5 mm
2	1	-	TBM-603	60 mm ball bearing Z stage +/-5 mm
3	1	-	GOHT-60A85BMS	60 mm Steel Black Chrome Goni, 85 mm CH, +/- 3.5 deg
4	1	-	GOHT-60A60BMS	60 mm Steel Black Chrome Goni, 60mm CH, +/-4.5 deg
5	1	-	KSP-606M	60 mm Black Chrome Rotatio stage, 360 deg coarse, +/-5deg fine
6	1	-	BSP-70120	Adapter plate, mounts to Inch or metric tables