






Product Presentation Bases

Optical Base







Catalog Code	Part Number / Product Name	Products image	Feature	Size, Material etc.
6032	OBC Optical Baseplates (With attached grips)		*Flatness: 0.2mm *With attached grips	300x450 - 450x600mm Tapped Holes on Top Surface: M4, M6, 1/4-20UNC Aluminum, Black Anodized
6019	OBA-LH Large Optical Rails		*X-shaped Cross Section *Rail Carriers (CAA-L) Available for System Construction *Custom Lengths and Scale Installation Available	Length: 500 - 2000mm Aluminum, Black Anodized
6020	OBT-LH/LUU Large Low-Profile Optical Rails		*Lightweight, Low axis, width 100mm *Custom Lengths and Scale Installation Available	Length: 500 - 2000mm Aluminum, Black Anodized
6021	CAA-L Carriers for Large Optical Rails		*Easy to be Slid on OBA-LH or OBT-LH/LUU and Secured at Arbitrary Position	Aluminum, Black Anodized
6022	OBA-SH Medium Optical Rails		*X-shaped cross section with fins extending radially in four directions. *Rail Carriers (CAA-LS) Available for System Construction *Custom Lengths and Scale Installation Available	Length: 300 - 2000mm Aluminum, Black Anodized
6023	OBT-SH/SUU Large Low-Profile Optical Rails		*Lightweight, Low axis, width 50mm *Custom Lengths and Scale Installation Available	Length: 300 - 2000mm Aluminum, Black Anodized
6024	CAA-LS/MB Carriers for Medium Optical Rails		*Carriers mount on a medium optical rail, slide freely and can be locked in place at any location	Aluminum, Black Anodized
6027	MB Magnet Base		*Switching the lever ON/OFF *Holding Force 49N - 490N	

Optical Base/ Spacers

Catalog Code	Part Number / Product Name	Products image	Feature	Size, Material etc.
6034 6035 6036	MSP-25/40 /60 Aluminum Spacers		*Allow Adjustment of Stage Height *Available in 1mm Increments (Except 25 x 25mm Size)	25x25, 40x40, 60x60mm Aluminum, Black Anodized
6037	SP/UP Adjustable Adapter Kits		*Lightweight, High rigidity *Posts are made to desired length in 1mm increments	Strut Length: $11 \leq A \leq 100$ (101mm or longer will be quoted separately) Plate material: Aluminum, Black Anodized Strut: Stainless steel
6058~ 6092	SP-*** /BSP Adapters Base Plate		*A Variety of Adapter Plates Available (see below the table)	Customized Shapes or Sizes Available Aluminum, Black Anodized


		Mounting hole pattern and screw shape of products to be mounted on upper surface													
		M16 Screw for fastening posts, post holders, etc.	20M2 Fasten 25x25 stages, etc.	31.75 6-32UNC Fasten 40x40 (For INCH) stages, etc.	32M3 Fasten 40x40 stages, etc.	40M3 Fasten 50x60 stages, etc.	50M4 Fasten 60x60 stages, etc.	P.C.D.50M4 Fasten 60 rotary stages, etc. (KSP-786 etc.)	50M6 Fasten 65x65 stages, etc.	50.8 1/4-20UNC 6Fasten 65x65 (For INCH) stages, etc.	70M4 Fasten 80x80 stages, etc.	90M4 Fasten 100x100 stages, etc.	100M4 Fasten 120x120 stages, etc.	M6 Internal thread Fasten minor holders, etc. in place of posts.	
Tapped hole pattern of products to be mounted on under surface	General Adapters	32M3 40x40 stage, etc. upper surface screw pattern	SP-132	SP-132		SP-132	SP-132	SP-132							
		P.C.D.32M3 40 rotary stage, etc. upper surface screw pattern	SP-129 SP-132	SP-129 SP-132		SP-129 SP-132	SP-129 SP-132	SP-132							
		50M4 60x60 stage, etc. upper surface screw pattern	SP-117 SP-133	SP-133		SP-133		SP-117 SP-133			SP-102-1	SP-102-2	SP-102-3	SP-122	
		P.C.D.50M4 60 rotary stage, etc. upper surface screw pattern	SP-129	SP-103 SP-129		SP-103 SP-129	SP-103 SP-129								
		31.75 6-32UNC 40x40 (For INCH) stage, etc. upper surface screw pattern		SP-401EE SP-401UU											
		P.D.C.31.75 6-32UNC 40 (INCH) rotary stage, etc. upper surface screw pattern		SP-401EE SP-401UU	SP-401UU	SP-401EE									
		50M6 65x65 stage, etc. upper surface screw pattern		SP-402EE	SP-402UU	SP-402EE			SP-402EE	SP-402UU					
		P.C.D.50M6 65/69 rotary stage, etc. upper surface screw pattern		SP-402EE	SP-402UU	SP-402EE			SP-402EE	SP-402UU					
		50.8 1/4-20UNC 65x65 (For INCH) stage, etc. upper surface screw pattern		SP-402EE	SP-402UU	SP-402EE			SP-402EE	SP-402UU					
	P.D.C.50.8 1/4-20UNC 65/69 rotary stage (For INCH), etc. upper surface screw pattern		SP-402EE	SP-402UU	SP-402EE			SP-402EE	SP-402UU						
	Special Adapters	MB-CB-CB	SP-130	SP-130		SP-130	SP-130	SP-130	SP-130						
		P.C.D.30M3	SP-129 SP-132	SP-123 SP-129 SP-132		SP-123 SP-129	SP-123 SP-129	SP-123 SP-129	SP-132						
		TAR-34*** TAR-38***D	SP-109-1 SP-109-2 SP-131			SP-131	SP-131	SP-131	SP-131						
		TAR-34403L TAR-34805L	SP-131			SP-131	SP-131	SP-131	SP-131						
		TAS-2060*	SP-127-1 SP-134	SP-127-1 SP-134		SP-134	SP-134	SP-134	SP-134						

Posts

Catalog Code	Part Number / Product Name	Products image	Feature	Size, Material etc.
6039	PST Pedestal Bases		*Connectable Design: M4 for Metric Type and 8-32 UNC for Inch Type *Easy to be Secured with Post Clamp	12.5 - 76.2mm Stainless steel, Finish: None
6040	PSC/PSCA/ PSCS/PSTM Pedestal Clamps		*Secure PST on Optical Breadboards or Tables *Available in Metric and Inch Types	Aluminum, Black Anodized (only PSCA) Stainless steel, Finish: None
6041	SP-T Adapter Plates for Pedestal		*For use with PST	Thickness: 5mm 25x25, 40x40mm Aluminum, Black Anodized
6042	PS-SP Pedestal Spacers		*For use with SP-T	Height: 0.5 - 4.3mm Stainless steel, Finish: None
6079	SPS Small Pedestal Base		*For use with $\phi 15$ Optics Holders *Low axis	Height: 12.5 - 35mm Stainless steel, Finish: None
6080	SPS-SP Small Pedestal Adapters		*Spacers for SPS	Thickness: 0.5 - 4mm Stainless steel, Finish: None

Pedestal Bases | PST


<Connection Example>




Holder Side	Supplied Screw	Post Stand Side
M6	CHS-M4-M6-9	M4
M4	SC-M4-18EE	M4
1/4"	CHS-1/4-8U-9	8-32UNC
8-32	SC-8U-18UU	8-32UNC

Small Pedestal Bases | SPS

<Connection Example>



<SPT/SPS Usage Example>



Posts / Poles

Catalog Code	Part Number / Product Name	Products image	Feature	Size, Material etc.
6052	RO/ROU /ROC Posts		*Stainless steel posts *External Thread on One End: M6 on RO-12/20, 1/4-20UNC on ROU, M4 on RO-6/ROC	φ6 - 20mm Length: 15 - 500mm Stainless steel, Finish: None
6083	RO-EE/UU Posts		*Connectable Design: M4 or M6 for Metric Type, 8-32UNC for Inch Type	φ12.7mm Length: 25.4 - 304.8mm Stainless steel, Finish: None
6053	PO Poles		*External M16P1 Thread on One End	φ12, 20mm Length: 100 - 500mm Stainless steel, Finish: None
6043	BRS High Stability Ball Plunger Post Holders		*High Stability Achieved by Two Lines and One Point Holding Design *Optional Base Available *Slide-Down Prevention Mechanism	φ12, 12.7, 20mm Height: 25.4 - 152.4mm Aluminum, Black Anodized
6044	RS Post Holders		*Value price *M4 or M16P1 External Mounting Thread at Bottom allows Mounting to Various Bases	φ6, 12, 20mm Height: 15 - 130mm φ6: Aluminum, Black Anodized φ12,20: Brass, Black matte
6081	RCA/RC Stand Base/Stand Clamp		*The base allows for posts to be positioned close together with the magnet of the clamps.	Aluminum, Black Anodized Stainless steel, Finish: None
6048	RSN Helicoidal Post Holders		*Height Adjustment via Rotation Ring *Approximately 5mm Travel per Revolution	φ12, 12.7, 20mm Aluminum, Black Anodized
6055 6054	CRHN/CCHN Cross Clamps		*Rotary and fixed type	φ6, 12, 12.7, 20mm Aluminum, Black Anodized

Vacuum (Bases)

High vacuum compatible products are now available as standard products.

Outgas from the material itself can be a problem in the ultra-high vacuum. Our products have selected the best material by examining the amount and type of outgas.


- **Stainless steel (SUS316L)**.....Strong corrosion resistance when exposed to high temperature baking, and less outgassing.
- **Aluminum (A6061)**.....Special surface treatment is used to reduce the degradation of the surface and outgassing.





The other parts also use the material that is suitable in the ultra-high vacuum.

■ Factor and product comparison by the degree of vacuum

Environment	Clean room	Low-medium vacuum	High vacuum	Ultra-high vacuum
Degree of vacuum	10 ⁵ Pa	10 ⁻¹ Pa	10 ⁻⁵ Pa	10 ⁻⁹ Pa
Exhaust system	Fan	Mechanical pump	Diffusion pump	Turbo-molecular pump Ionization pump, etc.
Major factor	Cleanliness	Residual gas	Volatile substance	Surface adsorption
General aluminum parts	○ Need Cleaning	△ Need vent holes	× Outgas from anodizing	× Baking generate contamination
Stainless steel parts	○ Need Cleaning	△ Need vent holes	△ Need vacuum grease	× Baking generate contamination
Vacuum compatible products	◎ Already cleaned	◎ Vent holes pre-designed	◎ Vacuum grease pre-applied	◎ Baking applicable

■ Cleaning



Catalog Code	Part Number / Product Name	Products image	Feature	Size, Material etc.
6108	VPST/VPSM Vacuum Pedestal Bases		*High Vacuum Compatible Post Stand *Low Outgassing *VPST: φ2 Vented body, VPSM: With screw holes for degassing	Height: 12.5 - 45mm SUS316L, Finish: None
6109	VPS-SP Vacuum Pedestal Spacer		*For use with VPST *Low Outgassing *Groove for Discharging Gas	SUS316L, Finish: None
6104	VPSCS Vacuum Pedestal Clamps		*High Vacuum Compatible Pedestal Clamp *Low Outgassing	SUS316L, Finish: None
6107	VRO Vacuum Posts		*High Vacuum Compatible Posts *φ2 Vented body	Length: 30 - 100mm SUS316L, Finish: None

Catalog Code	Part Number / Product Name	Products image	Feature	Size, Material etc.
6100	VRS Vacuum Post Holders		<ul style="list-style-type: none"> *For use with VRO *High Vacuum Compatible *$\phi 2$ Vented body 	Height: 30 - 100mm SUS316L, Finish: None
6106	VRCA Vacuum Stand Base/Stand Clamp		<ul style="list-style-type: none"> *For use with VRC *High Vacuum Compatible Post Holder Base *Low Outgassing *Groove for Discharging Gas 	Compatible: VRS-12, VRC-10 SUS316L, Finish: None
6105	VRC Vacuum Stand Clamps		<ul style="list-style-type: none"> *For use with VRC *High Vacuum Compatible Post Holder Base *Low Outgassing *With screw holes for degassing 	SUS316L, Finish: None
6101	VBSP Vacuum Base Plate		<ul style="list-style-type: none"> *High Vacuum Compatible Base Plate *With screw holes for degassing 	SUS316L, Finish: None
6102	VOBC Vacuum Optical Baseplates		<ul style="list-style-type: none"> *Low Outgassing compatible Optical Baseplates *With screw holes for degassing 	120x120mm - 200x300mm Aluminum (A6061) , Finish: None
6103	VBSC Vacuum Base Clamps		<ul style="list-style-type: none"> *Used for Fixing VOBC or VBSP *Fixable Thickness up to 10mm 	Secured with VPCB/VAGB Aluminum (A6061) , Finish: None
6111	VPCB/VAGB SDC® Clean bolt / Hexagon socket head cap screw		<ul style="list-style-type: none"> *Vacuum bolt with holes for degassing 	VPCB : SUS316L, SDC® Plasma surface hardening VAGB : SUS304, Finish: Silver plat
6110	VSC/VCHS Vacuum Dog Point Set Screws / Vacuum Thread Adapters		<ul style="list-style-type: none"> *For vacuum *Hexagon Socket Set Screw with Dog Point and M4/M6 Thread Adapter Available *Compatible with VPST and VRO 	SUS304, Finish: None