

Dear Valued Customers

**Notice for model change in Holders for
Laser Beam Expanders (BE-M)**

Thank you very much for your continued patronage.

Please be advised that new model of Holder for Laser Beam Expander (KLH-BE) has come into sale. New model's function is enhanced.
Accordingly we are terminating the production and sales of Holders for Laser Beam Expanders (BE-M).

In addition to the lens tilt adjustment function, centering mechanism is equipped with in the New model (KLH-BE).
Locking mechanism is also equipped with for keeping adjusted orientation.
Please refer to our web catalog for detailed specifications of new model, KLH-BE.

For part numbers of successor items, please refer to the following table.
We can also provide products of old models as custom made.
Please feel free to contact our sales division.

We appreciate your kind understanding and constant support.

Discontinued item (Old Part Number)	Replaced item (New Part number)	Compatible mounting Screw Size
BE-M22H	KLH-BE-M22H MHL-30BP※	M22 P0.75
BE-M34H	KLH-BE-M34H MHL-50BP※	M34 P1

※MHL-BP is base plate for MHL/KLH series. Post (RO-20-60) is not attached to the base plate. If Post Holder is used with, you can directly attach Rod (separately sold) to KLH-BE series Holders without using base plate.

Sincerely,