

Dec/2013

Dear Valued Customers

SIGMA KOKI CO., LTD.
Sales department,

Notice of the renewal of the motorized stage cable

Thank you very much for your continued business with Sigma Koki Co., Ltd.

Please be advised that we will start selling new motorized stage cable "D15RP cable" and "D15D15A cable" from January 2014 and accepting an order from December 2013.
For the specifications of a detail of the cable, please refer to our WEB catalog.

Along with this, sales of the "DMINIS cable", "DMINIH cable" and "DBCS cable" for the motorized stage will be the end of supply as of the end of February 2014.

For compatibility of the cable due to the sale end, please check the following correspondence table.

We apologize for any inconvenience this may cause.

Thank you very much again for your continued support.

[Notice]

New Products (starting sales from January/2014)

D15RP Cable

Model number	Price JPY [¥]	Connector (Controller side)	Connector (Stage side)	Length of cable [m]
D15RP-CA-2	6,500	17JE23150 Male Type	RP-17-13PA-12PC 12pin	2
D15RP-CA-3	7,000	17JE23150 Male Type	RP-17-13PA-12PC 12pin	3
D15RP-CA-5	9,000	17JE23150 Male Type	RP-17-13PA-12PC 12pin	5

D15D15 Cable

Model number	Price JPY [¥]	Connector (Controller side)	Connector (Stage side)	Length of cable [m]
D15D15A-CA-2	7,000	17JE23150 Male Type	17JE13150 Female Type	2
D15D15A-CA-3	9,000	17JE23150 Male Type	17JE13150 Female Type	3
D15D15A-CA-5	13,000	17JE23150 Male Type	17JE13150 Female Type	5

Discontinued

DMINIS Cable

DMINIS-CA-1
DMINIS-CA-2
DMINIS-CA-3
DMINIS-CA-4
DMINIS-CA-5

DMINIH Cable

DMINIH-CA-1
DMINIH-CA-2
DMINIH-CA-3
DMINIH-CA-4
DMINIH-CA-5

DBCS Cable

DBCS-15-2
DBCS-15-3
DBCS-15-4
DBCS-15-5

Correspondence table

Discontinued parts number		New Products parts number
DMINIS-CA-2	→	D15RP-CA-2
DMINIS-CA-3	→	D15RP-CA-3
DMINIS-CA-5	→	D15RP-CA-5
DMINIH-CA-2	→	D15RP-CA-2
DMINIH-CA-3	→	D15RP-CA-3
DMINIH-CA-5	→	D15RP-CA-5
DBCS-15-2	→	D15D15A-CA-2
DBCS-15-3	→	D15D15A-CA-3
DBCS-15-5	→	D15D15A-CA-5

If there are requests, we will be able to provide custom-made cable length of 1m or 4m.