

May 2015

Dear All Customers

Sigma Koki Co., Ltd.

### **Notification for specification changes in HPS/HDS series**

This is to notify that there have been some changes in specifications on Motorized Stages HPS/HDS series. In order to satisfy the requirements of CE marking, “High Performance Motorized Stages: HPS series” and “High Durability Motorized Rotation Stages: HDS series” will change the specifications. The changes will become effective from next production lot. Please refer to the tables below for the details.

Thank you for your continued support.

Part Number	Changes	Before	After
HPS series HPS60/80/120	Actuator	Uncovered	Covered
HDS-60YAW	Motor	TS3664N4 24 mm square (0.66A/Phase)	TS3664E10 24 mm square (0.75A/Phase)
HDS-80YAW	Dimension	<b>142</b> x 110.5 x 35 mm	<b>147</b> x 110.5 x 35 mm
	Motor	C9863-90215P 28 mm square (0.75A/Phase)	PK523HPB-C12 28 mm square (0.75A/Phase)
HDS-120YAW	Dimension	<b>164</b> x 161.9 x 35 mm	<b>169</b> x 161.9 x 35 mm
	Motor	C9863-90215P 28 mm square (0.75A/Phase)	PK523HPB-C12 28 mm square (0.75A/Phase)