

October, 2019

Dear customers,

SIGMAKOKI CO.,LTD.  
Sales Department

**Notice of Minor Change and Part Number**  
**Change**  
**(JD-100, JS-300, JB-400)**

Thank you very much for your constant patronage.

This is to notify that there will be changes in specifications on **Jog Dial (JD-100)**, **Joy Stick (JS-300)** and **Jog Operation Box (JB-400)** due to the end of main parts production by the manufacturer. The part numbers will be changed accordingly. Please visit our website for more details.

Please be advised that once remaining stock of current version of products is depleted, all quotation requests and orders for current version will be proceeded with for new version with new part number. Please visit our website for current stock status.

Please be notified that some controllers (see below for details) need version upgrade to be compatible with new version of above products.

Please feel free to contact us if you have any questions or concerns.

Thank you for your continued support.



Sincerely,

■Products with Minor Change

Product Name	Discontinued Part Number	Successor Part Number	Catalog Code
Jog Dial	JD-100	JD-101	W9085
Joy Stick	JS-300	JS-301	W9050
Jog Operation Box	JB-400	JB-401	W9050

■ Specification Changes

1. for Jog Dial, Joy Stick and Jog Operation Box:

Changes	JD-100 / JS-300 / JB-400	JD-101 / JS-301 / JB-401
Body Color	 <p>Cream</p>	 <p>White</p>

2. for Jog Dial only:

Changes	JD-100	JD-101
Speed Adjustment	Available for SHOT-30X only	Available for SHOT-30X, SHOT-702, PGC-04-U, HSC-103

3. for Jog Operation Box only:

Changes	JB-400	JB-401
Operable Axes	Axis 1/2	Axis 1/2/3/4

\* JS-301 does not have other specification changes except body color.

■ Regard the Controller Version Upgrade

JD-101, JS-301 and JB-401 are compatible with the following controller versions. If you are using controller with older version, please contact us for free version upgrade at our factory.

Controller Part Number	Catalog Code	Version
SHOT-302GS/304GS	W9047	Ver1.20~
SHOT-702	W9045	Ver1.05~
PGC-04-U	W9084	Ver1.02~
HSC-103	W9094	Ver1.01~