

Dear Customers,

SIGMA KOKI CO., LTD.  
Sales Department

## Notice of the Part Number Change

Thank you very much for your constant patronage.

Please be advised that the part numbers of the **Gimballed Beamsplitter Mounts (MHA-100S, MHA-130S)** and the **Dove Prism Holders (DBH-10)** will be changed to reflect improvements in the specifications.

The new parts include all the functionality of the old parts.

The following tables show the specification changes.

The current versions will no longer be sold once stocks have been depleted. The new versions will be launched in the middle of November 2016.


If you request/order the current model, we will quote the newest models depending on old stock quantities.

Current stock status can be seen on the web site.

Please do not hesitate to contact us if you have any questions.

Sincerely,  
OptoSigma Sales Team

Product name	Gimballed Beamsplitter Mounts	
	Current Model	New Model
Part Number	MHA-100S	MHA-100SA
Compatible Optics Thickness	4-14.5mm	4-15mm
C	20mm	21mm
Part Number	MHA-130S	MHA-130SA
Compatible Optics Thickness	7-17mm	7-18mm
C	23mm	24mm



Product name	Dove Prism Holders	
	Current Model	New Model
Part Number	DBH-10	DBHN-10
The scale plate	Fixed scale plate	Variable scale plate
A	φ59mm	φ60 mm
D	39mm	41mm
The lever	with	without

