

Dear Valued Customers,

Notice of end of production and successor products of MHF-50F (One Touch Mirror Mount)

Thank you very much for your constant patronage.

Please be announced that **MHF-50F (One Touch Mirror Mount)** will be discontinued at the end of February 2018.

Considering that its specification is overlapped with MHF-50.8F which can be used for both dia. 50 mm and 50.8 mm optics, we decided the discontinuation of MHF-50.

MHF-50.8F will be a replacement for MHF-50F in the beginning March, 2018.

When you want to put a post with 8-32 UNC screw on MHF-50.8F, please convert a screw size by AND-M6-8U.

Should you have any questions or concerns, please do not hesitate to contact us.

Thank you for your continued support.

Sincerely

Discontinued products	Successor products
MHF-50F	MHF-50.8F
<p>optics clamp (Push clamp) Switching mirror diameter by 90° rotation (68)</p> <p>Compatible optics $\phi 50.8$ (50)</p> <p>8-32UNC</p>	<p>Adjustment screw</p> <p>M6 P1</p> <p>2-$\phi 4.5$ through hole $\phi 8$ counterbore</p> <p>Mirror clamp (Push clamp) (Rotate the knob to switch to mirror size.)</p> <p>$\phi 50 / \phi 50.8$ (Mirror diameter)</p>