

June 2015

Dear All the Customers

Sigma Koki Co., Ltd.



### **Notification for specifications changes in MPH-12.7R**

Thank you very much for your constant patronage of our products.

This is to notify that there will be changes in specifications on polarizer holder MPH-12.7R/MPH-12.7R-N. Polarizer Holders MPH series are being updated to new models for the purpose of product improvement. MPH-12.7R/MPH-12.7R-N will be replaced by new model MPHN-12.7R/MPHN-12.7R-N in July 2015. For your information, the price will not be changed.

According to this change, old product MPH-12.7R/MPH-12.7R-N will be discontinued. If you would like to purchase of old model, please contact to International Sales Division.

Please refer to below of changes in specifications.

	Old Model	New Model
Part Number	MPH-12.7R/MPH-12.7R-N	MPHN-12.7R/MPHN-12.7R-N
Changes		
	Minimum Scale Resolution: 5°	Minimum Scale Resolution: 2.5°