

April, 2019

Dear Valued Customers,

SIGMAKOKI CO., LTD.
Sales Department

Notice of specification change of optical baseplate (OBC)

Thank you very much for your constant patronage.

Please be announced that the specification of **optical baseplate (OBC)** will be changed due to the soaring cost of materials.

The following specifications will be changed from the next production in order. There will be no change in the part number or price.

Please do not hesitate to contact us if you have any questions.
Thank you for your continued support.

Sincerely,

■ Optical Baseplate Catalogue Code: **W6032**

<Products>

Part Number		
Top surface tapping		
M4	M6	1/4-20UNC
OBC-3045	OBC-3045-M6	OBC-3045-1/4UU
OBC-4545	OBC-4545-M6	OBC-4545-1/4UU
OBC-4560	OBC-4560-M6	OBC-4560-1/4UU

<Description of change>

	Before	After
Surface treatment	Black anodized satin	Black anodized
Position of marking	Side	Top

<Position of marking image (reference) >

