

October 2015
Sigma Koki Co., Ltd.

Dear All the Customers

Notice of minor change of PGC-04

Thank you for your continued patronage.

Please be advised that there will be minor changes in specifications on **Pulse Generating Controllers PGC-04**.

Please refer to below of changes in specifications.

The changes will become effective from next production lot, scheduled in December 2015.

According to the change, old product PGC-04 will be discontinued.

	Old Model	New Model
Part number	PGC-04	PGC-04- U
Dimension (W x H x D)	195 x 140 x 60mm	180 x 140 x 60mm
Interface	RS232C / Ethernet	USB [*] / Ethernet

^{*}USB: USB port can be recognized as virtual serial port from external PC