

June 2015

Sigma Koki Co., Ltd

Dear Valued Customers

Notice of Discontinuation of 2 Mega-pixels USB Camera (SKDCE-3)

Thank you very much for your continued patronage.

Please be advised that we are terminating the sales of 2 Mega-pixels USB Camera (SKDCE-3) after our current stock is out, due to the high difficulty of gathering the parts.
For your reference, other 2 Mega-pixels USB Camera "[SK-TC202USB-AT](#)" is available.



Interface	USB2.0
Imaging Device	1/1.8 (2 Mega-pixels)
Chip Size	7.1×5.4mm
Pixels	1628×1236
Cell Size	4.4×4.4μm
Frame Rate	15ps
Lens Mount	C mount
Accessory	Viewer Software
	(USB camera driver Image loading software)
Operating System (OS)	Windows®XP、 Vista、 7 (32bit/64bit)
External Dimensions (excluding connector/protrusion)	51×51×46.1mm
Weight	0.19 kg

Contact our Sales Division for further information if required.

Sincerely,