

March 2016  
Sigma Koki Co., Ltd.  
Sales Department

Dear Valued Customers

**Notification for specification changes in TAMM series**

Thank you very much for your constant patronage of our products.

This is to notify that there have been some changes in specifications on Motorized Stages, TAMM series in order to satisfy the requirements of CE marking.

Please refer to the tables below for the details.

If you have any questions, please feel free to contact our sales division.

Thank you for your continued support.

Part Number	Changes	Old	New
TAMM40-10C/CR	Outer Dimension	<b><u>128.5</u></b> ×54.5×18※mm	<b><u>133.5</u></b> ×54.5×18※mm
	Motor	<b><u>C9863-90215P</u></b>	<b><u>PKP523HPB-C12</u></b>
TAMM60-15C/CR	Outer Dimension	<b><u>153.5</u></b> ×74.5×20※mm	<b><u>158.5</u></b> ×74.5×20※mm
	Motor	<b><u>C9863-90215P</u></b>	<b><u>PKP523HPB-C12</u></b>

※The motor height is 28mm