

March, 2017

Dear Customers,

Sigma Koki Co., Ltd.  
Sales Department

**Notice for the End of Production and  
Successor product of VNDA-90C1.5**

Thank you very much for your constant patronage.

Please be advised that we are terminating the sales of **Variable Reflective  
Aluminum ND Filter (VNDA-90C1.5)**.

As a successor product of VNDA-90C1.5, we have **Variable Reflective ND Filter,  
BK7 (VND-100)**.

For main differences, please refer to the following list.

We appreciate your kind understanding and constant support.

Sincerely,  
OptoSigma Sales Team

	Discontinued products	Successor products
Part Number	<b>VNDA-90C1.5</b>	<b>VND-100</b>
Transmittance	~98%	~92%
Size	Diameter φ90mm Inner diameter φ10mm Thickness 1.5mm	Diameter φ100mm Inner diameter φ20mm Thickness 3mm
Coating	Aluminum Dot pattern	Chrome Gradation
Parallelism	<3'	<1'
Price [JP Yen]	<b>40,000</b>	<b>80,000</b>