

March, 2017

Dear Customers,




SIGMA KOKI CO., LTD.
Sales Department

**Notice of end of production and
successor products of Microscope Systems (ZMM)**

Thank you very much for your constant patronage.

Please be advised that **Microscope Systems (ZMM-45**)** came to the end of sale and the new version, **Microscope Systems (ZMS-***)** have been launched as successor products.

Please do not hesitate to contact us if you have any questions.

Appearance			
New part number	ZMS-216	ZMS-316	ZMS-316H
Old part number	ZMM-45B1	ZMM-45T1	ZMM-45B2 ZMM-45T2
Type	Binocular / with stand	Trinocular / with stand	Trinocular / without stand
Total Magnification	7~45x(Zoom type)		
Working Distance(W.D.) [mm]	100		
Focusing Rack & Pinion Travel [mm]	47		
Illumination	not included	not included	not included
Attachment	Stand	Stand and C mount adapter	C mount adapter
Weight [kg]	4	4.3	2.9